



Alcohol free Disinfection wipes

For Surfaces and Medical Devices



30
secs

Effective against NOROVIRUS in **30 secs** Cleans and disinfects in one

PRODUCT SPECIFICATIONS

Sani-Cloth AF Universal

are disposable disinfection wipes for non-porous hard surfaces and non-invasive medical devices. Sani-Cloth AF Universal wipes can be used to pre-clean or decontaminate medical devices prior to sterilisation or high level disinfection.

	Test Method		Contact Time
Bacteria	EN 16615	<i>Pseudomonas aeruginosa</i>	60 seconds
	EN 16615	<i>Staphylococcus aureus</i>	60 seconds
	EN 16615	<i>Enterococcus hirae</i>	60 seconds
	EN 13727	<i>Pseudomonas aeruginosa</i>	30 seconds
	EN 13727	<i>Staphylococcus aureus</i>	30 seconds
	EN 13727	<i>Enterococcus hirae</i>	30 seconds
Enveloped Virus	EN 14476	HBV, HCV, HIV, H5N1, SARS, Corona	30 seconds
	EN 14476	Vaccinia Virus - Strain Elstree	30 seconds
Non-enveloped Virus	EN 14476	Norovirus	30 seconds
Yeasts	EN 16615	<i>Candida albicans</i>	60 seconds
	EN 13624	<i>Candida albicans</i>	30 seconds
TB	EN 14348	<i>Mycobacterium terrae</i>	2 minutes

- Multipurpose
Cleans and disinfects in one
- Broad spectrum of kill with low contact times
- Cost effective via product standardisation
- EN 16615**
Supported by extracted liquid testing
- Alcohol Free
PHMB Free
- Sani-Care**
Training Package*
*Terms & Conditions Apply



Class IIa
Medical Device

Products kill microbes but people prevent infections

BRO0005 Rev 3

200 Count Flow Wrap **Product code:** XP00390 **NHS code:** VJT592
Wipe size: 267mm x 210mm

40 Count Flow Wrap **Product code:** XP00391 **NHS code:** VJT640
Wipe size: 200mm x 200mm

100 Count Flow Wrap **Product code:** XP00363 **NHS code:** VJT521
Wipe size: 280mm x 200mm

200 Count Canister **Product code:** XP00366 **NHS code:** VJT514
Wipe size: 200mm x 200mm

225 Count Bucket **Product code:** XP00365 **NHS code:** VJT513
Wipe size: 270mm x 267mm



Products kill microbes but people prevent infections